

**New
From SDI!**

Astra[★]GUARD™

INFECTION PROTECTION

Bacterial/Viral Filter for Pulmonary Function Testing

AstraGuard – Infection Protection and Economy

The AstraGuard Bacterial/Viral Filter provides virtually perfect *infection protection* for expiratory and inspiratory spirometry testing. Its single-patient use greatly reduces costly and time-consuming disinfection procedures.

This easy-to-use disposable mouthpiece incorporates a highly efficient filter that virtually eliminates cross-contamination from patient to patient *and* from patient to health care professional.




The AstraGuard Bacterial/Viral Filter is designed for all Astra Spirometers.




AstraGuard Features

- 99.9% Bacterial and viral efficiency*
- Low resistance to flow
- Small, non-threatening size
- Comfortable oval mouthpiece
- Reduces risk of cross-contamination
- Less time spent disinfecting equipment



Cross-section of the filter shows placement of interior media. Edge of media is bonded to the exterior housing to form a tight seal.



The filter media is tailored to fit inside every mouthpiece, providing an effective guard against cross-contamination.

The AstraGuard Filter is available as a kit with The Klip, our noseclip.
Call your local distributor for details.

*As tested by Nelson Laboratories, September, 2006.

Your dealer is:

sdi

DIAGNOSTICS

The Spirometry Source

Ten Hampden Drive, Easton, MA 02375

www.sdidiagnostics.com