


# PRACTICE SAFE SPIROMETRY

## Infection Protection For Your Patient, Your Environment ...and You!

SDI Diagnostics is a leading manufacturer of spirometry supplies with over 20 years experience in PFT filters.

### *Pulmo*GUARD S™



### Introducing the PulmoGuard S a member of our new PulmoGuard family

- **Feel safer and save money**
- **Unique, pleated filter media**
- **Clear plastic for 100% visual inspection**
- **Recyclable**
- **Convenient packaging**
- **99.99% Bacterial and viral efficiency**
- **Meets ATS standards for resistance - .862cm H<sub>2</sub>O/L/Sec at 14L/Sec**
- **Dead space 73 ml.**
- **Infection protection for patients and health care professionals**
- **Exceeds al ATS/ERS standards for spirometry**



The Spirometry Source



